



CHOLESTECH LDX®
CALIBRATION VERIFICATION MATERIAL- RESULT SHEET

Lot No.: 0146

Issue Date: 2010-05-07

Expiration Date: 2011-05-26

Institution:	_____	Date Analyzed:	_____
Contact Person:	_____	Cholestech LDX s/n:	_____
Address:	_____	Cassette Lot No.:	_____
City/State/Zip:	_____	Operator:	_____
Telephone:	_____		
Fax:	_____		

TRIGLYCERIDE RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				137.8	111	164
Level 2				201.3	163	240
Level 3				306.6	248	365
Level 4				472.9	382	563

TOTAL CHOLESTEROL RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				137.1	114	160
Level 2				197.5	164	231
Level 3				307.1	255	359
Level 4				389.4	324	455

GLUCOSE RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				87.9	71	105
Level 2				209.4	169	250
Level 3				314.7	254	375
Level 4				423.7	343	500

HDL CHOLESTEROL RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				26.4	21	32
Level 2				36.6	28	45
Level 3				53.8	42	66
Level 4				69.8	54	85

Streck Part Number for Artwork Result Sheet: 350514-1

FAX results to Technical Service (858) 695-9964 or mail to:
 Biosite Incorporated, 9975 Summers Ridge Road, San Diego, CA 92121 USA
 Attn: Technical Service