

**CHOLESTECH LDX<sup>®</sup>**  
**CALIBRATION VERIFICATION MATERIAL- RESULT SHEET**  
 Issue Date: 2008-11-03 Expiration Date: 2009-11-12

Lot No.: 8317

Institution: _____	Date Analyzed: _____
Contact Person: _____	Cholestech LDX s/n: _____
Address: _____	Cassette Lot No.: _____
City/State/Zip: _____	Operator: _____
Telephone: _____	
Fax: _____	

**TRIGLYCERIDE RESULTS (mg/dL)**

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				136.6	110	163
Level 2				194.4	157	232
Level 3				288.5	233	344
Level 4				455.3	368	542

**TOTAL CHOLESTEROL RESULTS (mg/dL)**

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				143.2	119	167
Level 2				206.3	171	241
Level 3				323.3	269	378
Level 4				414.1	344	484

**GLUCOSE RESULTS (mg/dL)**

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				80.5	65	96
Level 2				205.5	166	245
Level 3				304.9	247	363
Level 4				419.9	339	500

**HDL CHOLESTEROL RESULTS (mg/dL)**

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				25.9	20	32
Level 2				38.4	30	47
Level 3				51.8	40	63
Level 4				69.8	54	85

Streck Part Number: 350514

FAX results to Cholestech Technical Service (510) 732-7227 or mail to:  
 Cholestech Corporation, 3347 Investment Blvd., Hayward, CA 94545-3808  
 Attn: Cholestech Technical Service