

**CHOLESTECH LDX®
CALIBRATION VERIFICATION MATERIAL- RESULT SHEET****Lot No.: 9063****Issue Date: 2009-02-12****Expiration Date: 2010-03-04**

Institution: _____

Contact Person: _____

Date Analyzed: _____

Address: _____

Cholestech LDX s/n: _____

City/State/Zip: _____

Cassette Lot No.: _____

Telephone: _____

Operator: _____

Fax: _____

TRIGLYCERIDE RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				131.2	106	156
Level 2				195.4	158	233
Level 3				281.7	228	336
Level 4				406.5	329	484

TOTAL CHOLESTEROL RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				135.8	113	159
Level 2				189.1	157	221
Level 3				286.5	238	335
Level 4				366.7	305	429

GLUCOSE RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				80.6	65	96
Level 2				199.5	161	238
Level 3				305.8	247	364
Level 4				391.0	316	466

HDL CHOLESTEROL RESULTS (mg/dL)

Sample ID	Results		Average	Target Mg/dL	Range	
	#1	#2			Low	High
Level 1				27.4	21	33
Level 2				37.2	29	45
Level 3				52.2	41	64
Level 4				63.6	49	78

Streck Part Number: **350514**

FAX results to Cholestech Technical Service (510) 732-7227 or mail to:
Cholestech Corporation, 3347 Investment Blvd., Hayward, CA 94545-3808
Attn: Cholestech Technical Service

PROPRIETARY INFORMATION:

Data herein is not to be disclosed without permission of Cholestech Corporation.