

PN: 402 14712

Rev. B

Date of last revision:  
2009/07/10

4.0" width  
3.0" length

Languages:

English

German

Spanish

Italian

French

Portuguese

Danish: Forventede værdier

Greek

Swedish: Förväntade värden

Front

**Cholestech** **REF** 14-446  
**LDX** **LOT** C17 7841

**hsCRP Controls**

EXPECTED VALUES / Erwartete Werte / Valores esperados / Valori previsti / Valeurs attendues / Valores esperados / Forventede værdier / Αναμενόμενες τιμές / Förväntade värden

hsCRP Controls <b>LOT</b> 0905-108	$\bar{x}$ (Mean)	$3 \times \sigma$ (3 Stdev)
L1 0905019	0.66	0.33-1.20
L2 0905020	3.81	1.81-7.02

Back

U.S.A.:  
Manufactured for  
Cholestech Corporation,  
3347 Investment Blvd.,  
Hayward, CA 94545 U.S.A