

PN: 402 14712

Rev. B

Date of last revision:

2009/07/10

4.0" width

3.0" length

Languages:

English

German

Spanish

Italian

French

Portuguese

Danish: Forventede værdier

Greek

Swedish: Förväntade värden

Front

Cholestech REF 14-446
LDX LOT C17 7959B

hsCRP Controls

EXPECTED VALUES / Erwartete Werte / Valores esperados / Valori previsti / Valeurs attendues / Valores esperados / Forventede værdier / Αναμενόμενες τιμές // Förväntade värden

hsCRP Controls <small>LOT</small> 0905-108	\bar{x} (Mean)	$3 \times \sigma$ (3 Stdev)
L1 0905019	0.61	0.30-1.10
L2 0905020	3.51	1.65-6.47

Back

U.S.A.:
Manufactured for
Cholestech Corporation,
3347 Investment Blvd.,
Hayward, CA 94545 U.S.A