



Lipid Control

L1 L2



Exp: 2009-11-12

L1 (Red Cap)

L2 (Black Cap)

LOT 8317

LOT 83170505

LOT 83170506










Constituent	Expected Range		Expected Range	
<b>For all Lipid Cassettes</b>				
Total Cholesterol	136-192 mg/dL	3.53-4.96 mmol/L	217-305 mg/dL	5.62-7.90 mmol/L
HDL Cholesterol	28-44 mg/dL	0.73-1.14 mmol/L	51-81 mg/dL	1.33-2.08 mmol/L
Triglycerides	100-148 mg/dL	1.13-1.67 mmol/L	213-314 mg/dL	2.40-3.54 mmol/L
Glucose	103-152 mg/dL	5.74-8.46 mmol/L	236-348 mg/dL	13.11-19.32 mmol/L

Sample Setting

Set the sample type to SAMPLE = "Whole B." in the Cholestech LDX configuration menu.



Cholestech Corporation  
3347 Investment Blvd.  
Hayward, CA 94545 USA

 • For In Vitro Diagnostic Use	 • Catalog Number	 • Level 1 and Level 2	 • Consult instructions for use	 • Temperature Limitation
 • Biological Risks	 • Lot Number	<b>Professional Use Only</b> • Professional Use Only	 • Manufacturer	 • Use By