



**Results in 5 minutes.**



**Doing Good.**



**Doing Well.**

**CHOLESTECH** 

3347 Investment Blvd.  
Hayward, CA 94545 U.S.A.  
Tel 800 733.0404  
Fax 510 732.7227  
[www.cholestech.com](http://www.cholestech.com)



**D o i n g**

**B o t h .**

Now you can help those most at risk in your community while generating goodwill and revenue for your facility. The Cholestech LDX System gives you a fast, simple way to identify those in your community who are most at risk for heart disease. More than ever before, people are aware of the potential dangers. As a result, cholesterol tests are in high demand.

Testing in your community, or providing these services to local employers, builds health system relationships and generates referrals to your qualified physicians.



Diabetes is one of the most prevalent and deadly chronic diseases in the world today. And proper glycemic control is a critical factor in preventing complications associated with the disease. But the fact is, in nearly 6,000 diabetic patients surveyed, two-thirds failed to achieve the ADA (American Diabetes Association) A1C target. A1C testing with the Cholestech GDX System can play a vital role in creating positive change within your community by helping those afflicted with diabetes maintain proper long-term glycemic control.

The immediate feedback of A1C results at the time of testing results in significant improvement of glycemic control, making your facility the vital link in maintaining and improving the health of your diabetic population.



Don't just do one.



A1C in 5 Minutes



Full Lipid Profile or ALT in 5 Minutes

Do both.

## Complete Lipid Profile or ALT in **5** minutes

- CLIA-waived
- Simple to use
- Lab-accurate results
- Uses a single drop of blood
- Small/portable design
- ALT (GPT) to assess the potential liver side effects of cholesterol-lowering, diabetes and other drug therapies
- Includes Framingham cardiac risk assessment feature

A lipid profile is a measurement of the levels of all the cholesterol in blood, both "good" and "bad," plus triglycerides.

### Results in 5 minutes.




#### Stick.

Collect blood sample from fingerstick and add to test cassette.



#### Click.

Insert cassette and press the Run button.



#### Done.

Results ready to discuss in just 5 minutes.



#### Available test cassettes:

|   |  |                              |                           |                                    |                               |                       |
|---|--|------------------------------|---------------------------|------------------------------------|-------------------------------|-----------------------|
| Lipid Profile<br>(TC, HDL, LDL, VLDL, Triglycerides, TC/HDL Ratio)<br>Cat. No. 10-989 | Lipid Profile <i>plus</i> Glucose<br>(TC, HDL, LDL, VLDL, Triglycerides, TC/HDL Ratio, Glucose)<br>Cat. No. 10-991 | ALT (GPT)<br>Cat. No. 11-772 | TC/HDL<br>Cat. No. 10-987 | TC/HDL/ Glucose<br>Cat. No. 10-990 | TC/Glucose<br>Cat. No. 10-988 | TC<br>Cat. No. 10-986 |
|---|--|------------------------------|---------------------------|------------------------------------|-------------------------------|-----------------------|

**Cholestech**  
**LDX**<sup>®</sup>

RESULTS. WHEN AND WHERE YOU NEED THEM.



# Cholestech GDX System. A1C in **5** minutes

- CLIA-waived
- Simple to use
- Lab-accurate results
- Uses a single drop of blood
- Small/portable design
- A1C gives the 2- to 3-month average of blood glucose levels as a measure of long-term glycemic control
- Meets NGSP certification requirements

## Results in 5 minutes.



### Stick.

Collect blood sample from fingerstick and add to sample tube.



### Twist.

Turn cartridge clockwise when you see flashing light and hear "beep." Pour solution. Repeat 3 times.



### Done.

Results ready to discuss in just 5 minutes.

**Cholestech**  
**GDX**<sup>TM</sup>  
RESULTS. IN TIME. EVERY TIME.



## Easy to Test Anywhere

The Cholestech LDX and GDX Systems are ideal in any venue where you want to provide lifesaving information for your people and ensure the health of your community...face-to-face in just 5 minutes.

- Health Screening Services
- Public Health Departments
- Government Agencies
- Military
- Pharmacies

The Cholestech LDX and GDX Systems can help the people in your community *and* help your facility. And *both* will thank you.

