

# The #1 Point-of-Care Lipid Testing System measures hs-CRP\*



**Cholestech**  
**LDX**

## BENEFITS AT A GLANCE

### AIDS IN PROVIDING PREDICTIVE VALUE

- CDC/AHA: assessment tool for moderate risk patients
- Adjunct to lipid testing

### Patients

- Face-to-face counseling with physician
- No lab visit or venous blood draw – saves time and money

### Practice

- Efficiencies for staff – eliminates costly call backs to lab and patient
- Increased revenue – reimbursed by Medicare in most states

### Physician

- Immediate intervention and counseling – better patient outcomes
- Fewer follow up visits – more free time

# The Cholestech LDX System brings hs-CRP testing to the point of care.



## hs-CRP TEST DETAILS

Methodology	Immunoassay	Assay range	Whole blood: 0.3 –10 mg/L Serum: 0.3 – 8 mg/L
Sample volume, whole blood	50µL		
Sample volume, serum	40µL	Assay duration, min.	Whole blood: 7 Serum: 4
Control volume	40µL		

## RELATIVE RISK CATEGORY AND AVERAGE hs-CRP LEVEL

Low : <1.0 mg/L

Average: 1.0 to 3.0 mg/L

High: >3.0 mg/L

## CHOLESTECH LDX Product Listing and Ordering Information

Cat. No.	Product
10-959	Cholestech LDX System <sup>1</sup>
12-807	hs-CRP Test Cassettes
12-918	hs-CRP Control Levels 1 and 2 (1 vial Level 1 & 1 vial Level 2, 3 mL each)
13-037	hs-CRP Starter Kit <sup>2</sup>
13-184	ROM Upgrade Kit
10-991	Lipid Profile•GLU
10-989	Lipid Profile
12-788	ALT•AST (GPT•GOT)
10-990	TC•HDL•GLU
10-987	TC•HDL
10-988	TC•GLU
10-986	TC (Total Cholesterol)
12-712	Multi-Analyte Control Levels 1 and 2 (1 vial Level 1 & 1 vial Level 2, 2 mL each)
12-713	Multi-Analyte Control Levels 1 and 2 (3 vials Level 1 & 3 vials Level 2, 2 mL each)

**CHOLESTECH**

<sup>1</sup>Includes Analyzer, Printer Cable, User Manual, Quick Reference Guide, Starter Pack, Optics Check Cassette, Printer Cable, User Manual, Procedure Manual, and Training Video. <sup>2</sup>100 Capillary Tubes, 100 Capillary Plungers, 50 Pipette Tips, 1 MiniPet Pipette (40µL), 1 MiniPet Pipette (50µL)

\*hs-CRP tests are for use in labs for moderate complexity testing and results in 7 minutes.

Cholestech LDX is a registered trademark of the Inverness Medical group of companies. Cholestech manufactures the LDX System.  
© 2008 Inverness Medical. All Rights reserved. MKT13685 Rev.B

**inverness medical**

Physician Diagnostics Group  
3347 Investment Blvd.  
Hayward, CA 94545-3808 USA  
Tel 877.441.7440  
Fax 510 732.7227  
www.cholestech.com